

	Okay as found	Corrective Action/Remarks
I. Defibrillator Unit		
Clean, no spills, and casing undamaged		
2. Cables/Connectors		
Inspect for damaged insulation, frayed/broken wires, or bent connector pins Connectors engage securely		
3. Sensors (pulse oximetry, NIBP cuff and hose, temperature sensors)		
Inspect for signs of damage or excessive wear Connectors engage securely		
1. Paddles		
Clean, not pitted, or damaged. Switches operate freely Cables and connectors free of damage and engage securely		
5. Supplies and Accessories		
a. Therapy pads in sealed pouches (2 sets, not expired) b. Defib gel or gel pads c. ECG monitoring electrodes d. Alcohol wipes e. Razors/scissors f. Recorder paper		
6. Batteries		
a. Fully charged battery installed in unit b. Fully charged spare battery available		
7. Ready For Use Indicator shows Ready - No flashing or Do Not Use symbol		
3. Operational Checks		
A. Power Up Sequence		
Battery charge and AC indicators illuminate with AC power Battery charge and Visual Alarm indicators briefly Illuminate Calif Test passed		
B. Hands Free Defibrillation (Test with only battery power)		
a. CHECK THERAPY ELECTRODES and PADS SHORT DETECTED messages display b. Charge time < 7 seconds c. DEFIB SHORT TEST passed at 30 Joules		
C. Pacer Test (Test with only battery power)		- "
a. Printer prints Pace markers every 25 mm at 60 ppm b. PACING: SHORT DETECTED message displays no error messages at 100 mA c. PACING: CHECK THERAPY ELECTRODES message displays		
D. Paddles		·······
a. APPLY PADDLES TO PATIENT message displays when paddles are connected to MFC b. Paddle switches functional (Recorder, Energy Select, Charge, Shock) c. DEFIB SHORT TEST passed at 30 Joules		
E. Reconnect Unit to AC Power		

Major problem(s) identified

Signature _

(OUT OF SERVICE)